



Required settings for Software Version 21.B

Starting from software version 21B, the scaling for the z-channels require a correct setting of piezo constants and High/Voltage Amplifier gains.

1 Changes in ini-file (HVA gains)

- Make yourself a high-level administrator by selecting and entering the password “anfatec”
- Open the used ini-file from the software:



- Search for the string “TestMic” in the ini file:

```
wave | erig | de | ay
[TestMic]
TestChannelComb
TestVolt0=0.5
TestVolt1=0.5
ADChannel=4
FormPosLeft=10
FormPosTop=10
DAScale0=9.8156
Gain0=7.0
```


- Add entries for the DA-channels 0 (this channel equals “Topo”, the entry is: “Gain0”), 2 and 3 (these channels equal x and y, the gain entries are “Gain2” and “Gain3”) as follows:

```
Gain0=15
Gain2=15
Gain3=15
GainzFast=15 // this is for the fast Z-channel in high-speed systems
```

These gains are the amplifications of the high-voltage amplifiers. In your case, these numbers should be 15.

2 Changes in the scanner file (piezo gains)

- Open the used scanner file from the software in Notepad with three steps:



Add the following lines:


```
XGain=0.39  
YGain=0.39  
ZGain=60  
ZFastGain=30
```

The numbers behind these entries are the piezo constants in units like “nm/V” for z-directions and “µm/V” for x- and y-direction. The volts are the volts of the HV amplifier.

Calculation example: A piezo stack with 150 V maximum voltage supply and 9 µm maximum movement range has a piezo gain of $9 \mu\text{m} / 150 \text{V} = 60 \text{ nm/V}$.

3 Check the settings

For the two possible z-channels, the piezo gain and HVA gain are displayed in the Z-tab of the scanner window:




For the X-Y-direction, the gains are used to display the scan range calibration coefficients in units of “nm/V” as typically provided by the piezo supplier:

Goto to the tab “Scan Range” I the scanner window.

	Range [µm]	X [µm/digit]	Y [µm/digit]
1	6	2.358E-8	2
2	10	2.62E-8	2
3	20	2.896E-8	2
4	30	3.183E-8	2
5	40	3.389E-8	2
6	50	3.447E-8	2

Enter the correct maximum voltage that the HV amplifier output can supply I the field below “Max. Range [V]”. For a -150 V to 150 V amplifier, this voltage is 300 V. For the 0 V to 150 V amplifier, this voltage is 150V.



Chose the tick “nm/V” instead of “µm/digit”, and the entries I the scan range table are recalculated in nm/V and displayed in this unit.

The vertical axis should be approximately the value that was entered as “XGain” in the scanner file.